

Product datasheet for **TA398667**

Goat Monoclonal Antibody [Clone ID: RMG03]

Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | RMG03 |
| Applications: | ELISA, WB |
| Recommended Dilution: | ELISA: 0.01 ug/mL – 0.5 ug/mL WB (WB): 0.1 ug/mL-0.5 ug/mL |
| Reactivity: | Rabbit |
| Host: | Goat |
| Isotype: | IgG |
| Clonality: | Monoclonal |
| Immunogen: | Rabbit IgG |
| Specificity: | This antibody reacts to Rabbit IgG. No cross reactivity with human IgG, rat IgG, or mouse IgG. |
| Formulation: | 50% Glycerol/PBS with 1% BSA and 0.09% sodium azide |
| Concentration: | 1 mg/mL - lot specific |
| Purification: | Protein A affinity purified from an animal origin-free culture supernatant |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |



[View online »](#)