

## Product datasheet for **TA398669**

### Goat Monoclonal Antibody [Clone ID: RMG06]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RMG06
Applications:	ELISA
Recommended Dilution:	ELISA: 0.05ug/mL – 1ug/mL
Reactivity:	Mouse
Host:	Goat
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Mouse IgG
Specificity:	This antibody reacts to the Fc region of Mouse IgG, including IgG1, IgG2a, and IgG2b. It shows very slight reaction with mouse IgG3. No cross reactivity with IgM, IgA, IgE, human IgG, rat IgG, or rabbit IgG.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein G affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)