

Product datasheet for **TA398670**

Ighg1 Goat Monoclonal Antibody [Clone ID: RMG07]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RMG07
Applications:	ELISA
Recommended Dilution:	ELISA: 0.05ug/mL – 1ug/mL
Reactivity:	Mouse
Host:	Goat
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Mouse IgG
Specificity:	This antibody reacts to Mouse IgG, including IgG1, IgG2a, IgG2b, and IgG3. No cross reactivity with IgM, IgA, IgE, human IgG, rat IgG, or rabbit IgG.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein G affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	immunoglobulin heavy constant gamma 1 (G1m marker)
Database Link:	Entrez Gene 16017 Mouse P01868
Synonyms:	IGHG1



[View online »](#)