

Product datasheet for **TA398676**

IGHA1 Rabbit Monoclonal Antibody [Clone ID: RM124]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM124
Applications:	ELISA, ICC, IHC
Recommended Dilution:	ELISA: 50ng/well – 200ng/well (for Capture) 0.05ug/mL – 0.2ug/mL (for Detection) Immunocytochemistry (ICC): 0.5ug/mL-2ug/mL Immunohistochemistry (IHC): 0.1ug/mL-1ug/mL
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Human IgA
Specificity:	This antibody reacts to human IgA1, and very slightly cross reacts with IgA2. No cross reactivity with human IgG, IgM, IgD, or IgE.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	immunoglobulin heavy constant alpha 1
Database Link:	Entrez Gene 3493 Human P01876



[View online »](#)