

## Product datasheet for **TA398678**

### **IgD (IGHD) Rabbit Monoclonal Antibody [Clone ID: RM123]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	RM123
Applications:	ELISA
Recommended Dilution:	ELISA: 25ng/well – 200ng/well (for Capture) 0.01ug/mL – 0.1ug/mL (for Detection)
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Human IgD
Specificity:	This antibody reacts to human IgD. No cross reactivity with human IgG, IgM, IgA, or IgE.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	immunoglobulin heavy constant delta
Database Link:	<a href="#">Entrez Gene 3495 Human P01880</a>



[View online »](#)