

## Product datasheet for **TA398679**

### IgE (IGHE) Rabbit Monoclonal Antibody [Clone ID: RM122]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM122
Applications:	ELISA
Recommended Dilution:	ELISA: 10ng/well – 100ng/well (for Capture) 0.01ug/mL – 0.1ug/mL (for Detection)
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Human IgE
Specificity:	This antibody reacts to human IgE. No cross reactivity with human IgG, IgM, IgD, or IgA.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	immunoglobulin heavy constant epsilon
Database Link:	<a href="#">Entrez Gene 3497 Human P01854</a>



[View online »](#)