

Product datasheet for **TA398684**

IGHG4 Rabbit Monoclonal Antibody [Clone ID: RM217]

Product data:

Product Type:	Primary Antibodies
Clone Name:	RM217
Applications:	ELISA
Recommended Dilution:	ELISA: 50ng/well – 200ng/well (for Capture) 0.05ug/mL – 0.2ug/mL (for Detection)
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Peptide corresponding to the Fc region of Human IgG4
Specificity:	This antibody reacts to the Fc region of human IgG4. No cross reactivity with human IgG1, IgG2, IgG3, IgM, IgA, IgD, IgE, mouse IgG, rat IgG, or goat IgG.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	immunoglobulin heavy constant gamma 4 (G4m marker)
Database Link:	Entrez Gene 3503 Human P01861



[View online »](#)