

## Product datasheet for **TA398685**

### Tubulin (TUBA1B) Rabbit Monoclonal Antibody [Clone ID: RM113]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM113
Applications:	ChIP, ICC, IHC, IP, WB
Recommended Dilution:	Western Blot (WB): 1:1000 dilution Immunoprecipitation (IP): 1:200 dilution Chromatin IP (ChIP): 1:200 dilution Immunocytochemistry (ICC): 1:200 dilution Immunohistochemistry (IHC): 1:200 dilution
Reactivity:	All Species
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A peptide corresponding to the C-terminus of $\alpha$ -Tubulin
Specificity:	This antibody reacts to $\alpha$ -Tubulin, including Tubulin alpha-1A chain and Tubulin alpha-1B chain.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	N/A - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tubulin alpha 1b
Database Link:	<a href="#">P68363</a>
Synonyms:	K-ALPHA-1



[View online »](#)