

## Product datasheet for **TA398688**

### Rabbit Monoclonal Antibody [Clone ID: RM146]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM146
Applications:	FC, ICC, IHC, IP, WB
Recommended Dilution:	Western Blot: 0.1 ug/mL – 0.5 ug/mL Immunocytochemistry 0.5 ug/mL – 2 ug/mL Flow Cytometry: 0.5 ug/mL – 2 ug/mL Immunohistochemistry 0.1 ug/mL-1 ug/mL ELISA: 0.01 ug/mL-0.5 ug/mL
Reactivity:	All Species
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Mixture of a peptide with 6xHis-Tag at the N-terminus and another peptide with 6xHis-Tag at the C-terminus.
Specificity:	This antibody reacts to recombinant proteins containing the 6xHis-Tag or 10xHis-Tag fused to either the amino or carboxy terminus. No cross reactivity with other endogenous protein in mammalian or bacteria cells.
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)