

## Product datasheet for **TA398690**

### Rabbit Monoclonal Antibody [Clone ID: RM105]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	RM105
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 0.02 ug/mL-0.5 ug/mL Immunohistochemistry (IHC): 0.5 ug/mL-2 ug/mL Western Blot (WB) : 0.1 ug/mL-1 ug/mL
Reactivity:	All Species
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	KLH-PEG with terminal methoxy group
Specificity:	This antibody reacts to the methoxy group of Polyethylene glycol (PEG). It is specific for methoxypoly(ethylene glycol).
Formulation:	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
Concentration:	1 mg/mL - lot specific
Purification:	Protein A affinity purified from an animal origin-free culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.



[View online »](#)