

## Product datasheet for **TA398729L**

### CCN2 Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Rat, Rabbit
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	CCN2 / CTGF
Formulation:	< 0.001 EU/μg, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144.86 kDa
Database Link:	<a href="#">P29279</a>



[View online »](#)