

## Product datasheet for **TA398775L**

### DLK1 Monoclonal Antibody

#### Product data:

|                         |  |
|-------------------------|--|
| Product Type:           | Primary Antibodies   |
| Applications:           | Animal Model, ELISA, FACS, FN, Kinetics  |
| Reactivity:             | Human, Mouse, Cynomolgus   |
| Isotype:                | IgG1   |
| Clonality:              | Monoclonal   |
| Specificity:            | DLK1   |
| Formulation:            | < 0.001 EU/μg, determined by LAL method.   |
| Concentration:          | 1 mg/ml  |
| Purification:           | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column. |
| Conjugation:            | Unconjugated   |
| Stability:              | Shelf life: one year from despatch.  |
| Predicted Protein Size: | 150 kDa  |
| Database Link:          | <a href="#">P80370</a>   |



[View online »](#)