

## Product datasheet for **TA398778**

### IL5RA Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Rat, Rabbit
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	IL-5Ra/ CD125
Formulation:	< 0.001EU/ $\mu$ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	146.06 kDa
Database Link:	<a href="#">Q01344</a>



[View online »](#)