

Product datasheet for **TA398780**

PCSK9 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	PCSK9
Formulation:	< 0.001EU/ μ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	146.24 kDa
Database Link:	Q8NBP7



[View online »](#)