

## Product datasheet for **TA398842L**

### FCRL5 Monoclonal Antibody

#### Product data:

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Applications:           | Animal Model, ELISA, FACS, FN, Kinetics   |
| Reactivity:             | Human, Cynomolgus   |
| Isotype:                | IgG1  |
| Clonality:              | Monoclonal  |
| Specificity:            | FcRH5 / IRTA2 / CD307e  |
| Formulation:            | < 0.001EU/μg,determined by LAL method.  |
| Concentration:          | 1mg/ml  |
| Purification:           | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Conjugation:            | Unconjugated  |
| Stability:              | Shelf life: one year from despatch.   |
| Predicted Protein Size: | 145 kDa   |
| Database Link:          | <a href="#">Q96RD9</a>  |

[View online »](#)