

## Product datasheet for **TA398865**

### **TNFRSF13C Monoclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	Animal Model, ELISA, FACS, FN, Kinetics
<b>Reactivity:</b>	Human, Cynomolgus
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Specificity:</b>	TNFRSF13C / BAFFR / CD268
<b>Formulation:</b>	< 0.001EU/μg,determined by LAL method.
<b>Concentration:</b>	1mg/ml
<b>Purification:</b>	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
<b>Conjugation:</b>	Unconjugated
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Predicted Protein Size:</b>	146.44 kDa
<b>Database Link:</b>	<a href="#">Q96RJ3</a>



[View online »](#)