

Product datasheet for **TA398879L**

PDCD1 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Cynomolgus
Isotype:	IgG4SP
Clonality:	Monoclonal
Specificity:	PDCD1 / PD-1 / CD279
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144.04 kDa
Database Link:	<u>Q15116</u>



[View online »](#)