

Product datasheet for **TA398888**

TFPI Monoclonal Antibody

Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | Animal Model, ELISA, FACS, FN, Kinetics |
| Reactivity: | Human, Macaque, Rabbit |
| Isotype: | IgG4SP |
| Clonality: | Monoclonal |
| Specificity: | TFPI |
| Formulation: | < 0.001EU/μg,determined by LAL method. |
| Concentration: | 1mg/ml |
| Purification: | Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column. |
| Conjugation: | Unconjugated |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 145 kDa |
| Database Link: | P10646 |



[View online »](#)