

Product datasheet for **TA398895**

CD22 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Rat, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Siglec-2 / CD22
Formulation:	< 0.001EU/ μ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	MMAE
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	151.97 kDa
Database Link:	P20273



[View online »](#)