

Product datasheet for **TA398914**

CD274 Monoclonal Antibody

Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | Animal Model, ELISA, FACS, FN, Kinetics |
| Reactivity: | Human, Mouse, Cynomolgus, Dog, Rat, Rabbit |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Specificity: | B7-H1 / PD-L1 / CD274 |
| Formulation: | < 0.001EU/ μ g,determined by LAL method. |
| Concentration: | 1mg/ml |
| Purification: | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Conjugation: | Unconjugated |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 150 kDa |
| Database Link: | Q9NZQ7 |



[View online »](#)