

## Product datasheet for **TA398927L**

### **TREM1 Monoclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	TREM1 / CD354
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	143.66 kDa
Database Link:	<a href="#">Q9NP99</a>

[View online »](#)