

Product datasheet for **TA398978L**

IL6 Monoclonal Antibody

Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | Animal Model, ELISA, FACS, FN, Kinetics |
| Reactivity: | Human, Macaque |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Specificity: | IL-6 / IFN γ |
| Formulation: | < 0.001 EU/ μ g, determined by LAL method. |
| Concentration: | 1 mg/ml |
| Purification: | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column. |
| Conjugation: | Unconjugated |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 145 kDa |
| Database Link: | P05231 |



[View online »](#)