

## Product datasheet for **TA398986L**

### **PDCD1 Monoclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Cynomolgus
Isotype:	IgG4SP
Clonality:	Monoclonal
Specificity:	PDCD1 / PD-1 / CD279
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in 0.1M Glycine,0.2M NaCl,(tris),pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145 kDa
Database Link:	<a href="#">Q15116</a>



[View online »](#)