

## Product datasheet for **TA398987**

### IL1RL2 Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	IL-1RL2 / IL-36R
Formulation:	< 0.001EU/ $\mu$ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145 kDa
Database Link:	<a href="#">Q9HB29</a>



[View online »](#)