

Product datasheet for **TA398991L**

ANGPT2 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Rat, Cynomolgus, Rabbit, Dog
Isotype:	IgG4SP
Clonality:	Monoclonal
Specificity:	ANGPT2
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.24 kDa
Database Link:	Q15123



[View online »](#)