

Product datasheet for **TA398993**

Integrin α 2 β 3 (ITGA2 & ITGB3) Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Integrin α 2 β 3 (ITGA2 & ITGB3)
Formulation:	< 0.001EU/ μ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in 0.1M Glycine,0.2M NaCl,(tris),pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.3 kDa
Database Link:	P17301 & P05106



[View online »](#)