

## Product datasheet for **TA398995**

### IL1B Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Macaque, Cynomolgus
Isotype:	IgG2SA
Clonality:	Monoclonal
Specificity:	IL-1b
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.12 kDa
Database Link:	<a href="#">P01584</a>



[View online »](#)