

Product datasheet for **TA399015L**

TIM-3 / HAVCR2 / CD366 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	TIM-3 / HAVCR2 / CD366
Formulation:	< 0.001 EU/μg, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145 kDa
Database Link:	Q8TDQ0-1

[View online »](#)