

Product datasheet for **TA399033L**

Integrin α V β 6 (ITGAV & ITGB6) Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Integrin α V β 6 (ITGAV & ITGB6)
Formulation:	< 0.001 EU/ μ g, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.12 kDa
Database Link:	P06756 & P18564



[View online »](#)