

Product datasheet for **TA399063L**

Integrin α 4 β 7 (ITGA4 & ITGB7) Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Integrin α 4 β 7 (ITGA4 & ITGB7)
Formulation:	< 0.001 EU/ μ g, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144.3 kDa
Database Link:	P13612 & P26010

[View online »](#)