

## Product datasheet for **TA399070**

### Integrin Alpha IIb Beta 3 Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Integrin Alpha IIb Beta 3
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	47.6 kDa
Database Link:	<a href="#">P08514-1 &amp; P05106-1</a>



[View online »](#)