

## Product datasheet for **TA399081L**

### ITGAL Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Macaque
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Integrin $\alpha$ L / ITGAL / CD11a
Formulation:	< 0.001 EU/ $\mu$ g, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	146.14 kDa
Database Link:	<a href="#">P20701</a>

[View online »](#)