

Product datasheet for **TA399128**

KLRC1 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Macaque, Cynomolgus, Baboon
Isotype:	IgG4
Clonality:	Monoclonal
Specificity:	NKG2A / CD94
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145 kDa
Database Link:	P26715



[View online »](#)