

Product datasheet for **TA399170L**

MSTN Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Rat, Chicken, Turkey, Cow
Isotype:	IgG4
Clonality:	Monoclonal
Specificity:	GDF8 / Myostatin
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	142.88 kDa
Database Link:	<u>Q14793</u>



[View online »](#)