

## Product datasheet for **TA399239L**

### SELL Monoclonal Antibody

#### Product data:

|                         |                                                                                                                                 |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | Primary Antibodies                                                                                                              |
| Applications:           | Animal Model, ELISA, FACS, FN, Kinetics                                                                                         |
| Reactivity:             | Human                                                                                                                           |
| Isotype:                | IgG2a                                                                                                                           |
| Clonality:              | Monoclonal                                                                                                                      |
| Specificity:            | L-selectin                                                                                                                      |
| Formulation:            | < 0.001EU/μg,determined by LAL method.                                                                                          |
| Concentration:          | 1mg/ml                                                                                                                          |
| Purification:           | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Conjugation:            | Unconjugated                                                                                                                    |
| Stability:              | Shelf life: one year from despatch.                                                                                             |
| Predicted Protein Size: | 145 kDa                                                                                                                         |
| Database Link:          | <a href="#">P14151</a>                                                                                                          |



[View online »](#)