

Product datasheet for **TA399266**

PDCD1 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Cynomolgus
Isotype:	IgG4SP
Clonality:	Monoclonal
Specificity:	PDCD1 / PD-1 / CD279
Formulation:	< 0.001EU/ μ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	146.54 kDa
Database Link:	Q15116



[View online »](#)