

## Product datasheet for **TA399274**

### EGFR Monoclonal Antibody

#### Product data:

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Applications:           | Animal Model, ELISA, FACS, FN, Kinetics   |
| Reactivity:             | Human, Cynomolgus   |
| Isotype:                | IgG1  |
| Clonality:              | Monoclonal  |
| Specificity:            | ERBB1 / EGFR / HER1   |
| Formulation:            | < 0.001EU/ $\mu$ g,determined by LAL method.  |
| Concentration:          | 1mg/ml  |
| Purification:           | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Conjugation:            | Unconjugated  |
| Stability:              | Shelf life: one year from despatch.   |
| Predicted Protein Size: | 145.28 kDa  |
| Database Link:          | <a href="#">P00533</a>  |



[View online »](#)