

Product datasheet for **TA399297L**

CCL20 Monoclonal Antibody

Product data:

| | |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | Primary Antibodies |
| Applications: | Animal Model, ELISA, FACS, FN, Kinetics |
| Reactivity: | Human, Cynomolgus, Macaque |
| Isotype: | IgG4 |
| Clonality: | Monoclonal |
| Specificity: | CCL20 |
| Formulation: | < 0.001EU/μg,determined by LAL method. |
| Concentration: | 1mg/ml |
| Purification: | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Conjugation: | Unconjugated |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 150 kDa |
| Database Link: | P78556 |



[View online »](#)