

Product datasheet for **TA399306L**

ITGAV Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG2SA
Clonality:	Monoclonal
Specificity:	Integrin α V / ITGAV / CD51
Formulation:	< 0.001 EU/ μ g, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144.58 kDa
Database Link:	P06756



[View online »](#)