

Product datasheet for **TA399307L**

APCS Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	Serum Amyloid P / SAP
Formulation:	< 0.001 EU/μg, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145 kDa
Database Link:	P02743



[View online »](#)