

Product datasheet for **TA399317**

IL1RL1 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Cynomolgus
Isotype:	IgG2SA
Clonality:	Monoclonal
Specificity:	IL-1RL1 / ST2 / IL-33R
Formulation:	< 0.001EU/ μ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.2 kDa
Database Link:	Q01638



[View online »](#)