

Product datasheet for **TA399320**

Rabies virus GP Monoclonal Antibody

Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | Animal Model, ELISA, FACS, FN, Kinetics |
| Reactivity: | Human |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Specificity: | Rabies virus GP |
| Formulation: | < 0.001 EU/μg, determined by LAL method. |
| Concentration: | 1 mg/ml |
| Purification: | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column. |
| Conjugation: | Unconjugated |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 144.44 kDa |



[View online »](#)