

Product datasheet for **TA399357**

Integrin $\alpha 5 \beta 1$ (ITGA5 & ITGB1) Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG4SP
Clonality:	Monoclonal
Specificity:	Integrin $\alpha 5 \beta 1$ (ITGA5 & ITGB1)
Formulation:	< 0.001EU/ μ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.4 kDa
Database Link:	P08648 & P05556



[View online »](#)