

## Product datasheet for **TA399364**

### **TNFRSF17 Monoclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	TNFRSF17 / BCMA / CD269
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	146.06 kDa
Database Link:	<a href="#">Q02223</a>



[View online »](#)