

## Product datasheet for **TA399379L**

### IL1RL1 Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Cynomolgus
Isotype:	IgG2SA
Clonality:	Monoclonal
Specificity:	IL-1RL1 / ST2 / IL-33R
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144.98 kDa
Database Link:	<a href="#">Q01638</a>



[View online »](#)