

Product datasheet for **TA399386L**

CD40LG Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Mouse, Macaque
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	TNFSF5 / CD40L / CD154
Formulation:	< 0.001EU/μg,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	145.66 kDa
Database Link:	P29965



[View online »](#)