

## Product datasheet for **TA399439L**

### IL22RA1 Monoclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Macaque, Cynomolgus
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	IL-22Ra
Formulation:	< 0.001 EU/μg, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	142.34 kDa
Database Link:	<a href="#">Q8N6P7</a>

[View online »](#)