

Product datasheet for **TA399562**

CSF1R Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human, Macaque, Cynomolgus
Isotype:	IgG4SP
Clonality:	Monoclonal
Specificity:	CSF1R / M-CSFR / CD115
Formulation:	< 0.001EU/ μ g,determined by LAL method.
Concentration:	1mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	147.18 kDa
Database Link:	P07333



[View online »](#)