

Product datasheet for **TA399566L**

CSF2 Monoclonal Antibody

Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | Animal Model, ELISA, FACS, FN, Kinetics |
| Reactivity: | Human, Mouse, Cynomolgus |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Specificity: | CSF2 / GM-CSF |
| Formulation: | < 0.001 EU/μg, determined by LAL method. |
| Concentration: | 1 mg/ml |
| Purification: | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column. |
| Conjugation: | Unconjugated |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 145 kDa |
| Database Link: | P04141 |

[View online »](#)