

Product datasheet for **TA399584L**

PF4 Monoclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	Animal Model, ELISA, FACS, FN, Kinetics
Reactivity:	Human
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	CXCL4 / PF4
Formulation:	< 0.001 EU/μg, determined by LAL method.
Concentration:	1 mg/ml
Purification:	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative. This antibody is purified through a protein A column.
Conjugation:	Unconjugated
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	144.5 kDa
Database Link:	P02776



[View online »](#)